Fusion XT OA Plus Unicompartmental Osteoarthritis Instructions

A ADJUST THE LATERAL (OUTSIDE) CONDYLE PAD

With brace on, walk 5 minutes to position brace properly and ensure that lateral condyle pad is applying lateral force to knee. If more lateral force is required to relieve symptoms turn upper and lower adjustment screws on side of hinge clockwise. This will move condyle pad inward (toward knee), applying more lateral force to knee.

- Upper and lower adjustment screws can be adjusted differently to accommodate various knee contours.
- For more comfort, different adjustments of upper and lower screws will distribute lateral force more evenly.
- Make sure large gel condyle pad is attached to the lateral (outside) condyle basket.

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Fitting Instructions for loss or gain in leg girth:

B ADJUSTING THE LATERAL CONDYLE PAD INWARD OR OUTWARD

If there is a loss or gain in leg girth (swelling) the pressure against the knee from the condyle pads may be too great or insufficient. The lateral condyle pad may be adjusted to accommodate a change in leg girth.

1 ADJUSTMENT FOR LOSS IN LEG GIRTH:

Use the tool provided to turn the upper and lower hinge adjustment screws clockwise. This will move the condyle pad inward, tighter against knee.

2 ADJUSTMENT FOR GAIN IN LEG GIRTH (SWELLING):

Use the tool provided to turn the upper and lower hinge adjustment screws counterclockwise. This will move the condyle pad outward, away from knee.

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